



# BODY BAG



江苏如翼医疗科技有限公司

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plaza, No19,

East Renmin Road,Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663

Whatsapp/Mobile:+ 86-18962287398;15962399164

Email:info@wanrooemed.com

www.wanrooemed.com    www.wanrooems.com



JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.



# Non-woven fabric body bag

## GA400A

Non-woven fabric is a new generation of environmental protection materials, with moisture-proof, breathable, flexible, light weight, no combustion-supporting, easy to decompose, non-toxic, non-irritating, colorful, low price, can be recycled and so on.

- 1. Sewing and processing.
- 2. Straight zipper or C-type zipper is available.
- 3. 4 children's webbing handles and 6 adults' webbing handles.
- 4. Waterproof and leakage proof.
- 5. #5 double end zipper.
- 6. Maximum loading capacity ≥120kg.
- 7. Suitable for camping, funerals and body handling.



No handle



with 4 handles



with 6 handles



PE coating inside

Product model	GA400A		
Model	Infant	Child	Adult
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)		
Material	Non-woven Fabric		
Size	70×50cm	170×80cm	215×80cm
Thickness	100gsm	110gsm	120gsm
Zipper	Straight Zipper Or C Type Zipper		
Handles	None	4	6
Maximum loading capacity	≥120kg		



# Non-Woven Fabric Body Bag(120gsm)

## GA400B

- 1. Blue body bags to improve the efficiency of forensic photographic research.
- 2. Type C zipper
- 3.4 external handles
- 4. Maximum bearing capacity ≥150kg
- 5. Special seam design on four sides to ensure no leakage of liquid and avoid infection risk.
- 6.#5 double end zipper
- 7. Zippered identification sunroof improves facial visibility.
- 8. Internal leakproof rubber layer is not easy to leak. Perfect for moving your body at any time.




Transparent window for easy identification

Product model	GA400B
Colors	DARK BLUE
Material	Non-woven Fabric
Size	225×90cm
Weight	120gsm
Zipper	C Type Zipper
Handles	4 Handles
Maximum loading capacity	≥150kg

# Non-Woven Fabric Body Bag(158gsm)

## GA400C

- 1. Dark blue is the international navy color. The body bag can be military.
- 2. Made of 128gPP and 30gPE, longer service life.
- This material is made of polyethylene and polypropylene, which has a longer service life. This bag remains elastic at low temperature.
- This allows it to be used in winter environments or morgues without fear of material "cracking".
- 4. Type C zipper
- 5. Waterproof and leakage proof
- 6. 6 external handles
- 7. PE coating is used inside to prevent infection.
- 8. #5 double end zipper
- 9. Maximum bearing capacity ≥180kg
- 10. Can be used in hospitals, funerals, war, disasters or emergencies.



PE coating inside

Product model	GA400C
Colors	DARK BLUE
Material	128gPP + 30gPE
Size	220×85cm
Weight	158gsm
Zipper	C Type Zipper
Handles	6 Handles
Maximum loading capacity	≥180kg



# PVC Body Bag(Economic)

## GA401

- 1. 0.2mm thickness of PVC material, edge using heat sealing processing to ensure 100% leakproof.
- 2. # 5 Double zipper head.
- 3. Straight zipper or C-type zipper is available.
- 4. There is no handle around the bag, and the weight of the product in a single room is ≥120kg.
- 5.Economical body bag, low price, worth buying.



Thickness: 0.2mm  
straight zipper



Thickness: 0.2mm  
C zipper



Thickness: 0.2mm  
C zipper

Product model	GA401 Economic		
Model	Baby	Child	Adult
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)		
Material	PVC		
Size	70×50cm	130×70cm	230×92cm
Thickness	0.2mm		
Zipper	Straight Zipper Or C Type Zipper		
Handles	None ( You can customize the handle if you want. )		
Maximum loading capacity	≥120kg		



PVC Body Bag(Medium Duty)

GA401

- 1. 0.2mm thickness of PVC material, edge using heat sealing processing to ensure 100% leakproof.
- 2. #5 Double zipper head.
- 3. C-type zipper design.
- 4. 4 children external webbing handles and 6 adults external webbing handles.
- 5. High bearing capacity of single product, bearing capacity ≥120kg.



Thickness: 0.2mm  
4/6 handles



Thickness: 0.2mm  
with underpad(Optional)

Product model	GA401 Medium Duty	
Model	Child	Adult
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)	
Material	PVC	
Size	130×70cm	230×92cm
Thickness	0.2mm	
Zipper	C Type Zipper(Other zippers are available)	
Handles	4 Handles	6 Handles
Maximum loading capacity	≥120kg	

PVC Body Bag(Heavy Duty)

GA401

- 1. 0.2mm thickness of PVC material, edge using heat sealing processing to ensure 100% leakproof.
- 2. No handle, bearing weight ≥ 120kg.
- 3. #5 double end nylon zipper.
- 4. Three-side zipper design.
- 5. Can be used below 5° (need to add refrigerant), can be refrigerated in the corpse freezer.



Thickness: 0.2/0.25mm  
three sides zipper

Product model	GA401 Heavy Duty
Size	230×92cm
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)
Material	PVC
Thickness	0.2mm
Zipper	Three sides straight zipper
Handles	None
Maximum loading capacity	≥120kg



# PVC Body Bag(Heavy Duty Three-Dimensional)

## GA4011

- 1. 0.25mm thickness of PVC material, side PP woven cloth sealing design.
- 2. #5 Double zipper head.
- 3. C-type zipper design.
- 4. Equipped with 10 sturdy external webbing handles, maximum loading capacity  $\geq 150\text{kg}$ .
- 5. Electronic heat seal structure complies with OSHA regulation 3130, leakage proof.
- 6. Transparent label bag with blue zipper for easy registration and transportation.
- 7. The cloverleaf double body bag creates a relatively large space for the recovery and transportation of the body. The placement of the body ensures respect for the deceased without causing wear and tear.



three-dimensional double layer



Product model	GA4011
Size	225×75×25cm
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request))
Material	PVC
Thickness	0.25mm
Zipper	C Type Zipper
Handles	10 Handles
Maximum loading capacity	$\geq 150\text{kg}$



# PVC Body Bag(Heavy Duty Three Layer)

## GA403

- 1.0.5mm thick PVC material with inner lining to prevent liquid leakage.
- 2. Internal lining 0.1 mm leakproof provides you additional protection in extreme cases.
- 3.4 external webbing handles for children and 6 external webbing handles for adults.
- 4.#5 Double zipper head.
- 5. C-type zipper design, easy to open.
- 6. Electronic heating seal in accordance with OSHA regulation 3130.
- 7. Maximum bearing weight ≥300kg.
- 8.Suitable for coroners, transport companies, funeral homes and medical examiners.



4 external webbing handles for children



0.5mm thick PVC material with inner lining to prevent liquid leakage.


Product model	GA403	
Model	Child	Adult
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)	
Material	PVC	
Size	150×100cm	230×110cm
Thickness	0.5mm	
Zipper	C Type Zipper	
Handles	4 Handles	6 Handles
Maximum loading capacity	≥300kg	



PVC Heavy Duty Body Bag (With Built-in Handles)

GA403A

- 1. 0.3mm/0.4mm thickness of PVC material, edge using heat sealing processing to ensure 100% leakproof.
- 2. #5 double zipper puller.
- 3. Straight zipper or C-type zipper is available.
- 4. Six heavy-duty handles with PVC hose inside, welded directly into the bag, and integrated into the material.
- 5. Maximum loading capacity ≥200 kg.



edge using heat sealing processing to ensure 100% leakproof.



Product model	GA403A		
Model	Child	Adult	Extra Large
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request))		
Process	Heat Seal Processing		
Size	150×100cm	220×100cm	250×120cm
Thickness	0.3mm	0.3mm	0.4mm
Zipper	Straight Zipper Or C Type Zipper		
Handles	4	6	6
Maximum loading capacity	≥200kg	≥200kg	≥250kg



PVC Heavy Duty Body Bag(View-able)

GA403B

- 1. 0.5mm thick PVC material, a total of three layers, containing lining to prevent liquid leakage.
- 2. The inner lining is made of 0.1mm leakproof transparent bag, which is not easy to leak.
- 3. PVC outer body bag and internal leakproof transparent bag provide protection for leakage.
- 4. Electronic heating seal in accordance with OSHA regulation 3130.
- 5. Unique internal optical transparent bag for viewing remains in the bag to minimize the risk of human exposure.
- 6. Double #5 type quality zipper.
- 7. Maximum loading capacity ≥300kg.
- 8. 8 Built-in handles that are welded directly inside the belt to fuse with the material for easy transportation during operation.



View-able, able to view the face of the dead without causing pollution

Product model	GA403B	
Model	Adult	Extra Large
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)	
Material	PVC	
Size	220×100cm	250×120cm
Thickness	0.5mm	
Zipper	C Type Zipper	
Handles	8 Handles	8 Handles
Maximum loading capacity	≥300kg	

PVC Heavy Duty Body Bag(With J Zipper)

GA403C

- 1. 0.5mm thick PVC material with lining to prevent liquid leakage.
- 2. Internal lining 0.1mm leakproof provides you additional protection in extreme cases.
- 3.10 webbing handles, secured by heavy rivets, increase load capacity, maximum load capacity ≥300 kg.
- 4. Unique J-shaped zipper that can clearly reveal the victim's identity just by pulling it open.
- 5.#5 Double end zipper.
- 6. Electronic heating seal in accordance with OSHA regulation 3130.
- 7. It can be used in morgues to store bodies, for funerals, or by official bodies such as police, military, emergency services.



J shape zipper, you can check the identity of the deceased just by pulling it a little

Product model	GA403C
Size	230×92cm
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)
Material	PVC
Thickness	0.5mm
Zipper	J Type Zipper
Handles	10 Handles
Maximum loading capacity	≥300kg




# PEVA Body Bag(With Absorbent Pad)

## GA4021

PEVA material is a very common environmentally friendly material with good softness and a little elasticity. PEVA material has low density and light weight. No seasonal, can resist ultra-low temperature, low temperature resistance 30 degrees easily through. Can be used for biodegradation; If it is not needed, whether it is discarded or burned, it is absolutely environmentally friendly and will not bring any harm to the environment.

- 1. The PEVA material is resistant to high and low temperatures and remains elastic at -25°C, which can be used in outdoor winter environments or morgue coolers.
- 2. Double-sided waterproof coating, high load-bearing strength, waterproof.
- 3. There are 4 (or 6-8) built-in handles that are welded directly into the body bag material. To bear weight of not less than 75-200 kg;
- 4. Rust-resistant # 5 double nylon zipper with choice of straight or U-shaped zipper design.
- 5. Electronic heating seal in accordance with OSHA regulation 3130.
- 6. Environmentally friendly, pollution-free, can be directly incinerated. Comply with OSHA and EPA combustion limits for chlorine gas release.



View-able, able to view the face of the dead without causing pollution



With windows, Velcro type(optional)

Product model	GA4021		
Model	Child	Adult	Extra Large
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)		
Process	Heat Seal Processing		
Size	150×100cm	220×100cm	250×120cm
Thickness	0.3mm	0.4mm	0.5mm
Zipper	Straight Zipper Or C Type Zipper		
Handles	4	6	8
Maximum loading capacity	≥180kg	≥250kg	≥250kg



PEVA Heavy Duty Body Bag(Three Layer)

GA4022

- 1. Adopt PEVA material with 0.3mm thickness and PP woven fabric on the back for reinforcement and suture to increase load bearing, suitable for cremation and/or burial.
- 2. #5 double zipper, C-type zipper design.
- 3. 10 external handles, sewn on both sides of the belt, can hold at least 250 kg.
- 4. Environmentally friendly, pollution-free, can be directly incinerated. Comply with OSHA and EPA combustion limits for chlorine gas release.
- 5. It remains elastic at -25°C, so it can be used in outdoor winter environments or morgue coolers without fear of "cracking".
- 6. You can use a marker pen to write messages on the plastic surface of the body bag. This often happens especially when collecting large numbers of bodies or in morgues and storing them in freezers.



10 braided strap handles  
PEVA+PP Woven fabric in the bottom



PP woven fabric on the back for  
reinforcement and suture to increase load bearing

Product model	GA4022
Size	245×115cm
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)
Material	PEVA+PP Woven fabric in the bottom
Thickness	0.3mm
Zipper	C Type Zipper
Handles	10 Handles
Maximum loading capacity	≥250kg



PEVA Body Bag(With Outside Carry Handles)

GA4024

- 1.PEVA material, not affected by blood and normal body fluids.
- 2. Suitable for cremation, chlorine-free, biodegradable, environmental protection.
- 3. Anti-rust #5 double nylon zipper, C-type zipper design.
- 4. Electronic heating seal in accordance with OSHA regulation 3130
- 5.4 or 6 sturdy webbing handles made from high quality high density woven material for easy handling.
- 6.Maximum loading capacity :≥120kg



4/6 high density woven Handles(can move)



with underpad(Optional)

Product model	GA4024		
Model	Baby	Child	Adult
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)		
Process	Heat Seal Processing		
Size	70×50cm	130×70cm	230×92cm
Thickness	0.2mm		
Zipper	Straight Zipper Or C Type Zipper		
Handles	None	4	6
Maximum loading capacity	≥120kg		




PEVA Body Bag(Communicable Disease Isolation)

GA404

- 1. Made of PEVA material, chlorine free, suitable for cremation and/or burial.
- 2. High quality #5 double side white zipper, full length C zipper.
- 3. All joints are electronically heated and sealed in accordance with OSHA regulation 3130.
- 4. Skylight with Velcro to improve facial visibility and prevent infection without contact with the deceased.
- 5. Special design of three sides can cover the zipper to prevent liquid leakage.
- 6. 6 External handles for adult body bags and 4 external handles for children body bags.
- 7. Bags for infectious disease carcasses (e.g. cholera and haemorrhagic fevers) and emergency situations to seal human carcasses prior to burial. Where possible, observe local customs regarding funeral rites.



unique optically clear window to allow safe viewing of remains in the bag.



Product model	GA404			
Model	Child		Adult	
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)			
Material	PEVA			
Size	130×70cm		220×100cm	
Thickness	0.2mm	0.3mm	0.2mm	0.3mm
Zipper	C Type Zipper			
Handles	4Handles	6Handles	4Handles	6Handles
Maximum loading capacity	≥120kg	≥200kg	≥120kg	≥200kg



# PEVA Chlorine Free Body Bag

## GA402

- 1. Made of PEVA material, chlorine free, suitable for cremation and/or burial.
- 2. #5 Double head white zipper with choice of straight zipper or C-type zipper design.
- 3. No handle design, no chlorine gas emissions from combustion, more suitable for cremation. Can be customized according to requirements.
- 4. The material remains elastic at minus 25 degrees, so it can be used in outdoor winter environments or morgue coolers without fear of "cracking".
- 5. Use a marker pen to write messages on the plastic surface of the body bag. This often happens especially when collecting large numbers of bodies or in morgues and storing them in freezers.
- 6. Environmentally friendly, pollution-free, can be directly incinerated. Comply with OSHA and EPA combustion limits for chlorine gas release.



View-able, able to view the face of the dead without causing pollution

Product model	GA402		
Model	Infant	Child	Adult
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)		
Process	Heat Seal Processing		
Size	70×50cm	130×70cm	230×92cm
Thickness	0.2/0.25mm		
Zipper	Straight Zipper Or C Type Zipper		
Handles	None (Custom handles are available)		
Maximum loading capacity	≥120/160kg		

# PEVA Chlorine Free Body Bag(Transparent)

## GA402

- 1. A translucent, self-sealing human remains.
- 2. Designed to replace traditional zippered white body bags as a more effective solution to prevent leakage and odors. This is especially true when the deceased begins to deteriorate or there is a risk of pathogens escaping.
- 3. Electronic heating seal in accordance with OSHA regulation 3130.
- 4. #5 double zipper, optional straight zipper or C-type zipper design.
- 5. Made of PEVA material, chlorine free, suitable for cremation and/or burial.
- 6. No handle design, no chlorine emissions from combustion, more suitable for cremation. Optional 4 or 6 handle designs.
- 7. Maximum loading weight ≥120kg.
- 8. It remains elastic at -25°C, so it can be used in outdoor winter environments or morgue coolers without fear of "cracking".
- 9. A marker can be used to write information on the plastic surface of the body bag, which often happens especially when collecting large numbers of bodies or in morgues and storing them in freezers.



A translucent, self-sealing human remains.

Product model	GA402 Transparent		
Model	Infant	Child	Adult
Colors	Transparency (Can available in other colors upon request)		
Process	Heat Seal Processing		
Size	70×50cm	130×70cm	230×92cm
Thickness	0.2/0.25mm		
Zipper	Straight Zipper Or C Type Zipper		
Handles	None (Custom handles are available)		
Maximum loading capacity	≥120/160kg		



# Woven Cloth Body Bag(PP Woven+PE)

## GA406B1

Woven fabrics are inexpensive, lightweight, environmentally friendly and recyclable. They are recognized as environmentally friendly products to protect the ecology of the earth. It has the advantages of good strength, breathable and waterproof, environmental protection, good flexibility, non-toxic and tasteless, and low price. If the material breaks down naturally outdoors, it has a life span of just 90 days. Indoors, it breaks down within eight years. And it is non-toxic, tasteless, no residue when burning, so it will not pollute the environment, so it can protect the environment.

- 1. Standard color: blue.
- 2. Do not sew or weld on the long side to avoid leakage and improve strength.
- 3. Using PP woven cloth and PE film material, the body bag is super anti-tear, which can protect the body well.
- 4. Made of PP woven fabric and PE film material, the body bag is super tear resistant and can protect the body well.
- 5. Connect 6 handles externally and fix them with heavy rivets to increase the weight.
- Maximum loading capacity ≥200 kg
- 6. #5 Double zipper, C-shaped zipper design.
- 7. Unlike vinyl bags, all materials remain flexible up to -36°C.



Connect 6 handles externally and fix them with heavy rivets to increase the weight.



High strength PP + PE special braided coating waterproof cloth

Product model	GA406B1
Size	230×92cm
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)
Material	150g PP + 30g PE
Weight	180g
Zipper	C Type Zipper
Handles	6 Handles
Maximum loading capacity	≥200kg



Woven Cloth Body Bag(PP Woven+PE,With 8 Handles)

GA406B2

- 1. Standard color: blue.
- 2. Do not sew or weld on the long side to avoid leakage and improve strength.
- 3. Made of PP woven fabric and PE film material, the body bag is super tear resistant and can protect the body well.
- 4. The interior is covered with a layer of PE material to prevent leakage and odor. PP material is degradable, beautiful and strong tolerance.
- 5. There are 8 handles in total. The sewing process is integrated with the body bag. Maximum load capacity ≥160kg.
- 6. #5 Double zipper, C-type zipper design.
- 7. Unlike vinyl bags, all materials remain flexible, up to -36°.



There are 8 handles in total.  
The sewing process is integrated with the body bag.



The interior is covered with a layer  
of PE material to prevent leakage and odor.

Product model	GA406B2
Size	230×92cm
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)
Material	150g PP + 30g PE
Weight	180g
Zipper	C Type Zipper
Handles	8 Handles
Maximum loading capacity	≥160kg



Woven Cloth Hunting Body Bag

GA406A

- 1. It is composed of woven cloth, with waterproof covering film, sealed together by sewing. This prevents blood from leaking into hard-to-clean areas and creating a bad smell.
- 2. The bottom and sides are 100% waterproof to prevent blood leakage into your vehicle. Can carry 150kg weight, very strong, long service life.
- 3. Each bag is equipped with a layer of non-woven waterproof film at the bottom to absorb blood and other moisture.
- 4. Keeps dust and bugs out, and makes it easier to hang meat to cool.
- 5. There are two ribbons, 50 cm long, easy to drag.
- 6. C-shaped zipper is more convenient for placing animal carcasses.
- 7. Suitable for single cab and double cab van. Custom bags are available on request for all vehicles and trailers.



The bottom and sides are 100% waterproof to prevent blood leakage into your vehicle.

Product model	GA406A
Size	160×100cm
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)
Material	Woven Fabric
Zipper	C Type Zipper
Handles	2 Handles
Maximum loading capacity	≥150kg

Pet Body Bag(For animals)

GA406PET

- 1. Using 150g PP+PE woven fabric material, it can not only in natural environmental conditions to reduce the environmentally harmless substances, but also has flame retardant, uv resistance and other functions.
- 2. The interior is covered with a layer of PE material to prevent leakage and odor. PP material is degradable, beautiful and strong tolerance.
- 3. Envelope zipper, double pull for tighter. Turn it on or off wherever needed.
- 4. It can be cremated safely. They do not contain PVC or chlorinated compounds, dioxins or furans.
- 5. Improve the euthanasia experience for animals, families, veterinarians and crematorium teams.
- 6. Completely leakproof and opaque, providing a quiet environment for your pet.



Special body bag for animals

Product model	GA406PET		
Model	Small	Medium	Large
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)		
Size	66×38cm	92×60cm	127×83cm
Weight	150g PP + PE		
Zipper	Straight Zipper Or C Type Zipper		
Handles	Not equipped with handle, can be customized.		



# Oxford Cloth Body Bag (With PVC Coating Inside)


## GA407

Oxford cloth is a kind of versatile and versatile fabric. It has the advantages of light and light texture, soft feel, good water resistance and good durability. Oxford cloth itself is strong and tough, wear resistance, fall resistance, pressure resistance, anti-scratch, easy to deal with casual collision and friction when traveling. Oxford cloth can still retain its original appearance after frequent and long-term use.

- 1. This body bag is suitable for funerals and transfer of bodies at any time.
- 2. It is made of Oxford cloth with PVC coating inside, with high bearing strength and waterproof.
- 3. There are 6 webbing handles around, with a load bearing of no less than 160-250kg.
- 4. All joints are sealed with electronic heating to contain body fluids and prevent leakage.
- 5. Double zipper head, linear or C-type zipper design can be selected.



All joints are sealed with electronic heating to contain body fluids and prevent leakage.



oxford cloth with PVC coating inside

Product model	GA407	
Model	210D	400D
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)	
Size	75×215cm	
Zipper	C Type Zipper	
Handles	6 Handles	
Maximum loading capacity	≥160kg	≥250kg



# Water Recovery Body Bag

## GA408

- 1. Made of highly visible nylon mesh material that is both visible in water and free to drain when removed from water. It is waterproof, light and tough.
- 2. 10 High strength handle design, easy to lift water and water.
- 3. Complete three-side zipper design, double-end #5 zipper, easy to restore in water.
- 4. High visibility International yellow with excellent underwater visibility.
- 5. floater grips for zipper pulls,allow divers to find and handle the zipper
- 6. Currently used by Coast Guard, Navy, police and Fire Department dive teams.
- 7. It is a very environmentally friendly material because it is easier to degrade.



Suitable for use in coast guard, US navy, police and fire department dive teams nationwide.



Unique Custom Floater Grips for zipper pulls to allow divers to find and handle closing the zipper.

Product model	GA408
Colors	Yellow / Orange
Size	230×92cm
Weight	180g
Zipper	Three Side Zipper
Handles	10 Handles
Maximum loading capacity	≥200kg



# PVC Vinyl-Coated Polyester Mesh Body Bag

## GA409

- 1. PVC sandwich cloth material, with good rain proof performance and good durability.
- 2. Connect 6 strap handles externally, secure with heavy rivets, and strengthen galvanized gasket. Increase bearing weight. Maximum loading capacity ≥200 kg.
- 3. The body bag is stitched by turning technique, which has the advantages of high density, high strength, no air permeability and no leakage.
- 4. #5 Double zipper, C-type zipper design, easy to put into the body.



PVC sandwich cloth material, with good rain proof performance and good durability.



Connect 6 strap handles externally, secure with heavy rivets, and strengthen galvanized gasket.

Product model	GA409
Size	230×92cm
Colors	White, Black, Green, Blue(Blue color to enhance post mortem forensic photography study without removing remains from the bag (Also available in black or white upon request)
Handles	6 Handles
Zipper	C Type Zipper
Material	PVC
Maximum loading capacity	≥200kg



Optional  
Handle



outside carry handles  
外置提手



built-in handles  
内置提手




xtra strong carry handles  
加固提手

ID tag



transparent ID bag  
透明ID袋



consisting of toe tag,  
exterior ID tag,  
personal effects ID tag.


body bag ID kit  
尸体袋标签

Underpad




underpad  
吸水垫

Zipper



straight zipper  
直拉链



C type zipper  
C型拉链



J type zipper  
J型拉链



Standard zipper  
标准拉链



U-shaped Double Zipper With Covering  
双边包裹型拉链



Metal ring zipper  
金属环拉链



Waterproof zipper  
防水拉链



#5 zipper head  
#5 拉链头



#8 zipper head  
#8 拉链头